


PZŁĄSER


CO₂ fractional laser machine

UNPRECEDENTED INNOVATION IN LASER TECHNOLOGY

With fractional laser resurfacing the laser beam is broken up or fractionated into many small micro beams which are separated so that when they strike the skin surface small areas of the skin between the beams are not hit by the laser and left intact. These small areas of untreated skin promote a much more rapid recovery and healing with less risk of complications. The small areas treated by the fractional micro beams, called micro treatment zones, cause sufficient laser injury to promote new collagen production and resultant facial skin rejuvenation.


Acne scars can be treated with the Fractional CO_2 laser in order to help resurface and recontour atrophic scars of the face, chest, and back. The fractional CO_2 laser has been used effectively in many countries to resurface acne scars and post-traumatic scars of the face and body.





FRACTIONAL Master

10600nm Fractional CO2 Laser



10600nm wavelength belongs to far-infrared radiation that has high absorption of water and generates heat, it absorbs too high water content as well as incises, evaporates and coagulates tissue with water.

WAVELENGTH(nm)

Adopts high power 60W RF type Laser generation imported from USA Access.


The emission energy is more stable, which can ensure the best effect. There is no need to add water to the RF tube, which is more convenient for maintenance, saves time and effort, and has a longer service life.



Comparison of RF tubes and ordinary glass tubes in the market:

Function comparison of RF-excited CO2, and normal Glass-tube CO2.

The Power




Glass tube laser:power is not stable.

Life Span



The RF-excited tube laser Co2 can be used at least 3 years, the glass tube laser Co2 can be used for one or two years.

Energy



Glass tube laser energy decrease dramatically from the beginning to the end.

Focal Spot Size





Glass Tube RF-excited Tube

Glass tube laser: spot size and spot density are different in size and shape



4 MODE WITH 9 APPLICATORS

for different target areas

Super Pulse Mode



"CUT"Applicator

Eliminating hypertrophy scars or depression and pigmentation.

3 Gyne mode



360° tightening the vagina and enhancing the elasticity and youthfulness.

2 Fractional Mode







Highly intelligent analysis of freckle and acne marks, stretch marks, skin rejuvenation, and scar repair.

4 Vulva Mode



Improving the texture and tone of the vulva, enhancing its lubrication and comfort.

MULTIPLE CUSTOM TREATMENT GRAPHICS ADJUSTABLE SHAPES, SIZES AND SPACING




















PZŁASER®

CO2 FRACTIONAL LASER MACHINE

Application:

- 1. Remove all kinds of wrinkles (brow veins, glabellas veins, neck veins, pregnancy veinsand so on)
- 2. Strengthen, tightening & whitening skin, face promotion, Skin rejuvenation, face lift.
- 3. Age spots and blemishes, acne and speckle treatment: blain acne, scabby acne, allergicacne, papilla acne, improving the appearance of liquidic skin and acne pit.
- 4. Treat all kinds of telangiectasia, blood vessels, blood silks.
- 5. Repair of sun damaged skin on the face, neck, shoulders and hands.
- Remove all kinds of neoplasms, wart, syringoma, sebaceous gland tumor, becker nevus, haemangioma, adipoma, steatoma, etc.
- 7. Vaginal tightening, vaginal rejuvenation.










FRACTIONAL LASER MACHINE

Easy&Simple GUI


One touch system for easy by intuitive button control





ADVANTAGES:

- 1. 60W RF laser tube.
- 2. USA Access METAL RF tube.
- 3. Fractional mode, normal, vulve, gyne 4 modes.
- 4. Korea imported seven joints light guiding arm.
- 5. Multiple scanning shapes.
- 6. Laser working mode: Continuous.
- 7. No pain, no cut, no side-effects.
- 8. No rebound! Apparently long-lasting result.
- 9. Particular designed intelligent laser control system ensures a precise treatment.
- 10. RF-excited Co2 laser energy shooting is an even steady and safe treatment process.
- 11. The best result for vaginal tightening vaginal rejuvenation.



SPECIFICATION

	Col fractional laser machine
	Co2 fractional laser machine
Nominal Wavelength	10600nm
Laser power	1-60W
Work mode	Fractional mode, normal, vulve, gyne
Pulse width	pulse mode: 1-100ms Super pulse mode: 0.1-10ms
Beam delivery	seven joint Articulated Arm
Power supply	AC 220V±10%; 50Hz±1Hz
Aiming beam	650nm
Smoke Exhausting System	Built-inairblow
Output spot area	20mm×20mm



























Tel: 0086 - 371 - 55677886 Fax: 0086 - 371 - 55677886 E-mail: info@pzlaser.com For more information: www.pzlaser.com

Headquarters: Building 13, Electronic Industrial Park, Hi-tech Development Zone, Zhengzhou, Henan, China