Before & After Plasma Treatment



PLASMA PEN IS THE WORLDWIDE #1 NITROGEN FIBROBLASTING DEVICE

Delivered with platinum standard, fully accredited training:


Our Plasma Pen converts electrical energy into an electrostatic energy which is transmitted to the applicator tip (probe) by impulse. Our device creates an electrical discharge between two electrodes separated by an 'insulating dielectric barrier'. A charge collects on the surface of

the barrier & discharges a flash of plasma via our 0.4mm probe – all within millionths of a second. Using a continuous direct energy source; we ionise nitrogen & oxygen in the air to create a plasma gas. This plasma gas is sustained by the unique fingertip and on-demand control of plasma our device offers.









Area's of use:



Why choose Plasma Pen?

Our Plasma Pen soft surgery sublimation:

- Treats, repairs, lifts, regenerates & rejuvenates almost any area of the skin
- Spectacularly improves & brightens skin tone, laxity & texture
- Significantly tightens loose skin & dramatically reduces & plumps out lines & wrinkles
- Delivers immediate & long lasting results typically 3 years without the need for fillers, injectables or invasive surgery



